



EasyPrep miRNA Extraction Kit

**For purification of total RNA, including miRNA
(and isolate miRNA from total RNA), from cells, tissues, and animal blood.**

Contents

Introduction.....	3
Kit Contents	4
Protocol	5

Introduction

EasyPrep miRNA Extraction Kit is designed for isolation and purification of miRNA, small interfering RNA (siRNA), small nuclear RNA (snRNA), and total RNA. Lysis buffer in the kit is optimized to exert excellent lysis ability and isolation sensitivity. The kit utilizes a silica-based system to enrich small RNA from various sample sources, especially for small RNA of <200 nt. This kit could be applied for RNA isolation from various samples (cell, animal tissue, plant tissue, serum, plasma). Each single spin column could handle 30-50 mg animal tissue (for RNA enriched tissue like liver, should be less than 30 mg), 100 mg plant tissue or 1×10^7 cells. RNA of no contamination of DNA and protein could be obtained within 1 hour, and used in Northern Blot, Dot Blot, Poly A screening, in vitro translation, RNase protection analysis and also molecular cloning.

Important Notes

1. Change gloves frequently. Bacteria on the skin can result in RNase contamination.
2. Use RNase-Free plastic and tips to avoid cross contamination.
3. RNA can be protected in Buffer TZ. But RNA must be stored or applied in RNase-Free plastic or glassware. To wipe off RNase, the glassware can be heated at 150°C for 4 hours, and plastic can be dipped in 0.5 M NaOH for 10 min, washed by RNase-Free ddH₂O thoroughly, and sterilized.
4. Use RNase-Free ddH₂O to prepare solution. (Add DEPC to 0.1% final concentration in ddH₂O. Shake solution to mix, and leave overnight at room temperature, autoclave for 15 min).
5. Add 48 ml and 18 ml ethanol (96-100%) to Buffer TW and Buffer TPD respectively before use.

Kit Contents

Contents	DPT-BE01 (50 preps)
Buffer TZ	60 ml
Buffer TW ^a	15 ml
Buffer TPD ^b	12 ml
RNase-free ddH ₂ O	15 ml
Tmi Columns	50
TR Columns	50
Collection Tubes(1.5 ml)	50
Collection Tubes (2 ml)	2x50

- a. Prepare TW by adding 48 ml of ethanol (96-100%).
- b. Prepare TPD by adding 18 ml of ethanol (96-100%).

Storage

EasyPrep miRNA Extraction Kit should be stored at room temperature (15–25°C).

Buffer TZ should be stored at 2-8°C protected from light.

Protocol

A. Purification of miRNA from tissue and cell.

*Apply this protocol when there is high demand of miRNA purity, e.g. miRNA chip and miRNA clone research.

1. Preparation of samples
 - a. Tissue: Grind tissue in liquid nitrogen. Add 1ml Buffer TZ for every 30-50 mg animal tissue (or 100 mg plant tissue), homogenize minced using a tissue homogenizer. Sample volume should not be over 1/10 of Buffer TZ.
 - b. Monolayer cell: Add 1 ml Buffer TZ per unit (10 cm²). Pipet to mix and ensure that no cell clumps are visible.

Note: Addition volume of Buffer TZ depends on monolayer area, not cell number. If Buffer TZ is not enough, RNA obtained will be contaminated with DNA.

- c. Cell suspension: Centrifuge for 5 min at 2,100 rpm (400 × g), carefully remove all supernatant by aspiration and disrupt the cells by adding 1 ml Buffer TZ, vortex or pipet to mix. Don't wash cells before adding Buffer TZ, otherwise mRNA will be degraded.
2. Place the tube containing the homogenate on the benchtop at room temperature (15-25°C) for 5 min, to separate nucleic acids and protein.
3. Optional: Centrifuge the lysate at 12,000 rpm (~13,400 × g) for 5 min at 4°C to remove any particulate material. Then transfer supernatant to a new tube.

Note: When preparing samples with high content of fat, proteins, polysaccharides, or extracellular material (e.g., muscle, fat tissue, or tuberous plant material), an additional centrifugation may be required to remove insoluble material from the samples. The precipitation after centrifugation contains cell outer membrane, polysaccharide, high molecular mass DNA, RNA is in the upper supernatant.

4. Add 200 µl chloroform to the supernatant and cap it securely. Shake the tube vigorously for 15 s, put the tube containing the homogenate on the bench top at room temperature (15-25°C) for 5 min.
 5. Centrifuge for 15 min at 12,000 rpm (~13,400 × g) at 4°C. After centrifugation, the sample separates into 3 phases: an upper, colorless, aqueous phase containing RNA; a white interphase; and a lower organic phase. The volume of aqueous phase is around 50% of Buffer TZ added, transfer the aqueous phase to a new tube.
 6. Add 0.43 volume of Ethanol (96-100%) (e.g. add 215 µl Ethanol (96-100%) to 500 µl transferred liquid) and mix thoroughly by pipetting up and down several times (A precipitate may form after addition of ethanol, but this will not affect the procedure). Transfer the obtained liquid, include any precipitate that may have formed into a Tmi Column, centrifuge at 12,000 rpm (~13,400 × g) for 30

EASYPREP MIRNA EXTRACTION KIT

- s at room temperature (15-25°C). If the volume exceeds 700 µl, centrifuge successive aliquots in the same column. Discard the Tmi Column after centrifugation, and keep the flow-through.
7. Add 0.75 volume of Ethanol (96-100%) (e.g. add 525 µl Ethanol (96-100%)) to 700 µl flow-through) and mix thoroughly by pipetting up and down several times (A precipitate may form after addition of ethanol, but this will not affect the procedure). Transfer the obtained liquid, include any precipitate that may have formed into a TR Column, centrifuge at 12,000 rpm (~13,400 × g) for 30 s at room temperature (15-25°C). If the volume exceeds 700 µl, centrifuge successive aliquots in the same column. Discard the flow-through, and keep the TR Column.
 8. Add 500 µl Buffer TPD (ensure that 18 ml of ethanol has been added) to the TR Column. Close the lid gently and incubate 2 min at room temperature (15-25°C). Then centrifuge for 30 s at 12,000 rpm (~13,400 × g) to wash the column. Discard the flow-through.
 9. Add 500 µl Buffer TW (ensure that 48 ml of ethanol has been added) to the TR Column. Close the lid gently and incubate 2 min at room temperature (15-25°C). Then centrifuge for 30 s at 12,000 rpm (~13,400 × g) to wash the column. Discard the flow-through.
 10. Repeat step 9.
 11. Place the TR Column into a new 2 ml collection tube, centrifuge at 12,000 rpm (~13,400 × g) for 1 min, and discard the flow-through.

Note: Perform this step to eliminate any possible carryover of Buffer TW. After centrifugation place the TR Spin Column at clean bench for a while, to completely dry the column membrane. Residual Buffer TW will have negative influence on following RT experiment.

12. Transfer the TR Column into a new 1.5 ml RNase-Free Centrifuge Tube, add 15-30 µl RNase-Free ddH₂O directly onto the TR Column membrane and incubate 2 min at room temperature (15-25°C). Close the lid gently and centrifuge for 2 min at 12,000 rpm (~13,400 × g) to elute the RNA. Note: The volume of elution buffer should not be less than 15 µl. Small elution volume may have a negative effect on RNA yield. Store RNA solution at -70°C. To obtain a higher total RNA concentration, please repeat step 12.

B. Purification of total RNA from tissue and cells (Obtained total RNA includes small RNA like miRNA).

*Use protocol when there is not high requirement of miRNA purity, e.g. miRNA RT-PCR, and miRNA Northern blot research.

1. Preparation of samples

- a. Tissue: Grind tissue in liquid nitrogen. Add 1ml Buffer TZ for every 30-50 mg animal tissue (or 100 mg plant tissue), homogenize minced using a tissue homogenizer. Sample volume should not be over 1/10 of Buffer TZ.
- b. Monolayer cell: Add 1 ml Buffer TZ per unit (10 cm²). Pipet to mix and ensure that no cell clumps are visible.

Note: Addition volume of Buffer TZ depends on monolayer area, not cell number. If Buffer TZ is not enough, RNA obtained will be contaminated with DNA.

- c. Cell suspension: Centrifuge for 5 min at 2,100 rpm (400 × g), carefully remove all supernatant by aspiration and disrupt the cells by adding 1 ml Buffer TZ, vortex or pipet to mix. Don't wash cells before adding Buffer TZ, otherwise mRNA will be degraded.

2. Place the tube containing the homogenate on the benchtop at room temperature (15-25°C) for 5 min, to separate nucleic acids and protein.
3. Optional: Centrifuge the lysate at 12,000 rpm (~13,400 × g) for 5 min at 4°C to remove any particulate material. Then transfer supernatant to a new tube.

Note: When preparing samples with high content of fat, proteins, polysaccharides, or extracellular material (e.g., muscle, fat tissue, or tuberous plant material), an additional centrifugation may be required to remove insoluble material from the samples. The precipitation after centrifugation contains cell outer membrane, polysaccharide, high molecular mass DNA, RNA is in the upper supernatant.

4. Add 200 µl chloroform to the supernatant and cap it securely. Shake the tube vigorously for 15 s, put the tube containing the homogenate on the bench top at room temperature (15-25°C) for 5 min.
5. Centrifuge for 15 min at 12,000 rpm (~13,400 × g) at 4°C. After centrifugation, the sample separates into 3 phases: an upper, colorless, aqueous phase containing RNA; a white interphase; and a lower organic phase. The volume of aqueous phase is around 50% of Buffer TZ added, transfer the aqueous phase to a new tube.
6. Add 1.5 volume of Ethanol (96-100%) (e.g. add 750 µl Ethanol (96-100%) to 500 µl transferred liquid) and mix thoroughly by pipetting up and down several times (A precipitate may form after addition of ethanol, but this will not affect the procedure).

Transfer the obtained liquid, include any precipitate that may have formed into a Tmi Column, centrifuge at 12,000 rpm (~13,400 × g) for 30 s at room temperature (15-25°C). If the volume exceeds

EASYPREP MIRNA EXTRACTION KIT

- 700 μ l, centrifuge successive aliquots in the same column. Discard the flow-through, and keep the Tmi Column after centrifugation.
7. Add 500 μ l Buffer TPD (ensure that 18 ml ethanol has been added) to the Tmi Column. Close the lid gently and incubate 2 min at room temperature (15-25°C). Then centrifuge for 30 s at 12,000 rpm (~13,400 \times g) to wash the column. Discard the flow-through.
 8. Add 500 μ l Buffer TW (ensure that 48 ml of ethanol has been added) to the Tmi Column. Close the lid gently and incubate 2 min at room temperature (15-25°C). Then centrifuge for 30 s at 12,000 rpm (~13,400 \times g) to wash the column. Discard the flow-through.
 9. Repeat step 8.
 10. Place the Tmi Column into a new 2 ml collection tube, centrifuge at 12,000 rpm (~13,400 \times g) for 1 min, and discard the flow-through.

Note: Perform this step to eliminate any possible carryover of Buffer TW. After centrifugation place the Tmi Spin Column at clean bench for a while, to completely dry the column membrane. Residual Buffer TW will have negative influence on following RT experiment.

11. Transfer the Spin Column miRspin to a new 1.5 ml RNase-Free Centrifuge Tube, add 30-100 μ l RNase-Free ddH₂O directly onto the Tmi column membrane and incubate 2 min at room temperature (15-25°C). Close the lid gently and centrifuge for 2 min at 12,000 rpm (~13,400 \times g) to elute the RNA.

Note: The volume of elution buffer should not be less than 30 μ L. Small elution volume may have a negative effect on RNA yield. Store RNA solution at -70°C. To obtain a higher total RNA concentration, please repeat step 11.

C. Purification of miRNA from whole blood, serum or plasma.

1. Preparation of samples Add equal volume of Buffer TZ to whole blood, serum or plasma, vortex 30 s to mix thoroughly.

Note: The starting volume of sample should be at least 200 μ l, otherwise RNA yield will be low.

2. Place the tube containing the homogenate on the benchtop at room temperature (15-25°C) for 5 min, to separate nucleic acids and protein.
3. Centrifuge the lysate at 12,000 rpm (\sim 13,400 \times g) for 10 min at room temperature (15-25°C) to remove any particulate material. Then transfer supernatant to a new RNase-Free centrifuge tube.
4. Add 200 μ l chloroform to the supernatant and cap it securely. Shake the tube vigorously for 15 s, incubate the tube containing the homogenate on the bench top at room temperature (15-25°C) for 5 min.
5. Centrifuge for 15 min at 12,000 rpm (\sim 13,400 \times g) at 4°C. After centrifugation, the sample separates into 3 phases: an upper, colorless, aqueous phase containing RNA; a white interphase; and a lower yellow organic phase. Transfer the aqueous phase to a new tube.
6. Add 1/3 volume of Ethanol (96-100%) (e.g. add 100 μ l Ethanol (96-100%) to 300 μ l transferred liquid) and mix thoroughly by pipetting up and down several times (A precipitate may form after addition of ethanol, but this will not affect the procedure). Transfer the obtained liquid, include any precipitate that may have formed into a Tmi Column, incubate at room temperature (15-25°C) for 2 min, then centrifuge at 12,000 rpm (\sim 13,400 \times g) for 30 s. Discard the Tmi Column after centrifugation, and keep the flow-through.
7. Add 2/3 volume of Ethanol (96-100%) (e.g. add 200 μ l Ethanol (96-100%) to 300 μ l flow-through) and mix thoroughly (A precipitate may form after addition of ethanol, but this will not affect the procedure). Transfer the obtained liquid, include any precipitate that may have formed into a TR Column, incubate at room temperature (15-25°C) for 2 min, centrifuge at 12,000 rpm (\sim 13,400 \times g) for 30 s. Discard the flow-through, and keep the Tmi Column.
8. Add 500 μ l Buffer TPD (ensure that 18 ml of ethanol has been added) to the TR Column. Close the lid gently and incubate 2 min at room temperature (15-25°C). Then centrifuge for 30 s at 12,000 rpm (\sim 13,400 \times g) to wash the column. Discard the flow-through.
9. Add 500 μ l Buffer TW (ensure that 48 ml of ethanol has been added) to the TR Column. Close the lid gently and incubate 2 min at room temperature (15-25°C). Then centrifuge for 30 s at 12,000 rpm (\sim 13,400 \times g) to wash the column. Discard the flow-through.
10. Repeat step 9.
11. Place the TR Column into a 2 ml collection tube, centrifuge at 12,000 rpm (\sim 13,400 \times g) for 1 min, and discard the flow-through.

Note: Perform this step to eliminate any possible carryover of Buffer TW. After centrifugation place the TR Column at clean bench for a while, to completely dry the column membrane. Residual Buffer TW will have negative influence on following RT experiment.

12. Transfer the TR column to a new 1.5 ml RNase-Free Centrifuge Tube, add 15-30 μ l RNase-Free ddH₂O directly onto the TR column membrane and incubate 2 min at room temperature (15-25°C). Close the lid gently and centrifuge for 2 min at 12,000 rpm (\sim 13,400 \times g) to elute the RNA.

Note: The volume of elution buffer should not be less than 15 μ L. Small elution volume may have a negative effect on RNA yield. Please store RNA solution at -70°C. To obtain a higher total RNA concentration, please repeat step 12; or increase sample volume and Buffer MZ and chloroform volume in proportion.

The product is for research only, not for diagnostic and clinical use.