



TOOLS Bacterial and Fungal DNA Extraction Kit

Contents

Introduction.....	2
Kit Contents	2
Workfolw.....	3
Protocol	4

Introduction

TOOLS Bacterial and Fungal DNA Extraction Kit is a buffer-based system for bacterial and fungal DNA extraction and purification. Simple centrifugation procedures enable the complete removal of contaminants and enzyme inhibitors. It provides a fast, simple, and cost-effective method of purifying DNA that is suitable for downstream applications, including PCR, southern blot, genomic DNA libraries, and sequencing.

Kit Contents

Contents	TX-BFD01 (50 preps)
Buffer BDA	30 ml
Buffer BDB	12 ml
Buffer BDC	300 μ l
Binding Gel	50 tubes

Materials not provided


1-Bromo-3-chloropropane (BCP; CAS Number: 109-70-6).

TOOLS Proteinase K (Cat. No. RTT-BD03) or other compatible reagents and lysozyme.

Storage

TOOLS Bacterial and Fungal DNA Extraction Kit can be stored at room temperature for up to 24 months.

Workfolw



Protocol

Sample preparation

1. Place 1 mL of bacterial cells into a 2-mL microcentrifuge tube.
2. Centrifuge for 5 min at 1,000 rpm to pellet the cells and remove as much supernatant as possible.
3. Add 500 μ L of Buffer BDA to the tube and vortex for 60 s to resuspend the pellet thoroughly.
4. Add 5 μ L of proteinase K (20 mg/mL) and 5 μ L of lysozyme (100 mg/mL) to the tube, vortex for 60 s, and incubate for 1 h at 55 °C. Vortex the tube for 20 s to mix thoroughly.

DNA extraction

1. Add 5 μ L of proteinase K and 5 μ L of lysozyme to the tube. Vortex for 60 s and incubate at 55 °C for 30 min.
2. Add 200 μ L of BCP and 200 μ L of Buffer BDB to the tube and mix the sample by inverting the tube three times. Transfer the solution to the binding gel tube (centrifuge at 13,000 $\times g$, 30 s before use).
3. Invert the binding gel tube three times (do not vortex) and centrifuge at 12,000–16,000 $\times g$ for 5 min.
4. Transfer the supernatant to a new microcentrifuge tube and add 5 μ L of Buffer BDC and 550 μ L of isopropanol to the tube.
5. Invert the tube three times and incubate at 37 °C for 10 min.
6. Centrifuge at 14,000 rpm for 5 min. Remove the supernatant.
7. Add 600 μ L of 70% ethanol. Centrifuge at 14,000 rpm for 5 min. Remove the supernatant.
8. Air dry the pellet and rehydrate it with 20–50 μ L of TE buffer or ddH₂O (adjust buffer volume according to pellet size).

This product is for research only. Not for diagnostic or clinical use.